

# Medium FixSepar<sup>®</sup> ECO

Biosepar is pleased to introduce our newly-developed formaldehyde-free fixer and carrier solution for use in parasitology. **Medium FixSepar<sup>®</sup> ECO** is a perfectly integrated addition to the **Biosepar ParasiTrap<sup>®</sup> Diagnostic System**.

This new development is based on a ground-breaking concept, demonstrating Biosepar's innovative strength.

The new solution has the following characteristics:

- not harmful to the environment
- non-toxic
- completely safe for work and transport

The European regulations demand the preferential use of reagents which pose no or only limited danger to living organisms or the environment. Biosepar has succeeded in being the first to offer this product on the international market.

Usage:

**Medium FixSepar<sup>®</sup> ECO** can be used as a drop-in replacement for Biosepar Standard Medium A, i.e. 3,5 ml per test (please refer to the manual and to “Tips for Optimal Use of the Biosepar ParasiTrap<sup>®</sup> Diagnostic System”) **Medium FixSepar<sup>®</sup> ECO** delivers the same high quality results that could be obtained with formaldehyde-containing Medium A.

The **ParasiTrap<sup>®</sup> Diagnostic System** offers the following advantages:

- time-saving work method
- fully standardized
- reproducibility
- high sensitivity and hit rate for all parasites
- optimized quality control and assurance before, during and after test preparation
- no danger of infection
- closed system resulting in prescribed workplace-values
- compatible with all standard laboratory equipment
- fulfils all diagnostic and hygienic recommendations (WHO, CDC, DPG, BGA etc.)
- virtually non-polluting

Should you have any questions, please get contact us right away. We are looking forward to your call. Experience success and high-quality diagnostic results with **ParasiTrap<sup>®</sup>** and **Medium FixSepar<sup>®</sup> ECO**.

**Yours Biosepar Team**

Head office  
Am Industriepark 3, D-84453 Mühldorf am Inn  
Export office Munich  
Otto-Heilmann-Str. 2, P.O.Box 1344, D-82027 Grünwald/Munich  
Telefon: +49 (0) 8631 167 48 55, Fax: +49 (0) 8631 167 48 57  
E-mail: mail@biosepar.de  
Internet: www.biosepar.de

