


biosepar *ANIMAL line*







**THE NEW DEVELOPED PARASITOLOGICAL
DIAGNOSTIC SYSTEM FOR THE VETERINARY MEDICINE**

FAECAL DIAGNOSTICS



SMALL ANIMALS

ParasiTrap® Flotation PLUS I SA

Order-No.: 7486 FLOT I SA		(10, 50, 100 pieces/box)	Complete system	
				
Processing tube I	Processing tube II	FixSepar® Flotation PLUS Medium	FixSepar® Flotation PLUS Medium	
(filled with 5,5 ml Medium A FixSepar® Flotation solution)	(with integrated filter system)	(in brown glass with graduated dropper)	(incl. the sufficient quantity of media)	
				
Sample spoon		Screw cap		

Contents of packagings

10 units:

- 10 pcs. Processing tube I, filled with 5,5 ml FixSepar® Flotation PLUS Medium
- 10 pcs. Processing tube II, with integrated filter-system
- 10 pcs. Sample spoon
- 10 pcs. Screw cap
- 1 pc. 30 ml FixSepar® Flotation PLUS Medium in brown glass with graduated dropper

50 units:









- 50 pcs. Processing tube I, filled with 5,5 ml FixSepar® Flotation PLUS Medium
- 50 pcs. Processing tube II, with integrated filter-system
- 50 pcs. Sample spoon
- 50 pcs. Screw cap
- 1 pc. 125 ml FixSepar® Flotation PLUS Medium
- 1 pc. 30 ml FixSepar® Flotation PLUS Medium in brown glass with graduated dropper

100 units:

- 100 pcs. Processing tube I, filled with 5,5 ml FixSepar® Flotation PLUS Medium
- 100 pcs. Processing tube II, with integrated filter-system
- 100 pcs. Sample spoon
- 100 pcs. Screw cap
- 1 pc. 300 ml FixSepar® Flotation PLUS Medium
- 1 pc. 30 ml FixSepar® Flotation PLUS Medium in brown glass with graduated dropper

SMALL ANIMALS

ParasiTrap® Premium CombiPack Comfort System FixSepar® ECO

Order-No.: 7487 ECO	(10, 50, 100 pieces/box)	Premium Combi System	
			
Processing tube I (filled with 3,5 ml Medium A FixSepar® ECO solution)	Processing tube II (with integrated filter system)	Brown glass (with graduated dropper)	Combi Medium (incl. the sufficient quantity of media)
			
Sample spoon	One-way pipette	Cotton wool stick	Screw cap

Contents of packagings

10 units:

- 10 pcs. Processing tube I, filled with 3,5 ml Medium A FixSepar® ECO solution
- 10 pcs. Processing tube II, with integrated filter-system
- 10 pcs. Sample spoon
- 10 pcs. Screw cap
- 10 pcs. One-way pipette 1 ml
- 10 pcs. Cotton wool stick
- 1 pc. 30 ml Combi Medium in brown glass
- 1 pc. Graduated dropper

50 units:








- 50 pcs. Processing tube I, filled with 3,5 ml Medium A FixSepar® ECO solution
- 50 pcs. Processing tube II, with integrated filter-system
- 50 pcs. Sample spoon
- 50 pcs. Screw cap
- 50 pcs. One-way pipette 1 ml
- 50 pcs. Cotton wool stick
- 1 pc. 75 ml Combi Medium in brown glass
- 1 pc. Graduated dropper

100 units:

- 100 pcs. Processing tube I, filled with 3,5 ml Medium A FixSepar® ECO solution
- 100 pcs. Processing tube II, with integrated filter-system
- 100 pcs. Sample spoon
- 100 pcs. Screw cap
- 100 pcs. One-way pipette 1 ml
- 100 pcs. Cotton wool stick
- 1 pc. 150 ml Combi Medium in brown glass
- 1 pc. Graduated dropper

SMALL ANIMALS

ParasiTrap® Premium Combipack Comfort System AF

Order-No.: 7488 AF		(10, 50, 100 pieces/box)		Premium Combi System	
					
Processing tube I (filled with 3,5 ml Medium AF solution)		Processing tube II (with integrated filter system)		Brown glass (with graduated dropper)	
					
Sample spoon		One-way pipette		Cotton wool stick	
					
				Screw cap	

Contents of packagings

10 units:

- 10 pcs. Processing tube I, filled with 3,5 ml Medium AF solution
- 10 pcs. Processing tube II, with integrated filter-system
- 10 pcs. Sample spoon
- 10 pcs. Screw cap
- 10 pcs. One-way pipette 1 ml
- 10 pcs. Cotton wool stick
- 1 pc. 30 ml Combi Medium in brown glass
- 1 pc. Graduated dropper

50 units:









- 50 pcs. Processing tube I, filled with 3,5 ml Medium AF solution
- 50 pcs. Processing tube II, with integrated filter-system
- 50 pcs. Sample spoon
- 50 pcs. Screw cap
- 50 pcs. One-way pipette 1 ml
- 50 pcs. Cotton wool stick
- 1 pc. 75 ml Combi Medium in brown glass
- 1 pc. Graduated dropper

100 units:

- 100 pcs. Processing tube I, filled with 3,5 ml Medium AF solution
- 100 pcs. Processing tube II, with integrated filter-system
- 100 pcs. Sample spoon
- 100 pcs. Screw cap
- 100 pcs. One-way pipette 1 ml
- 100 pcs. Cotton wool stick
- 1 pc. 150 ml Combi Medium in brown glass
- 1 pc. Graduated dropper

SMALL ANIMALS

ParasiTrap[®] Premium Combipack Comfort System SAF

Order-No.: 7489 SAF	(10, 50, 100 pieces/box)	Premium Combi System	
			
Processing tube I <small>(filled with 3,5 ml Medium A SAF solution)</small>	Processing tube II <small>(with integrated filter system)</small>	Brown glass <small>(with graduated dropper)</small>	Combi Medium <small>(incl. the sufficient quantity of media)</small>
			
Sample spoon	One-way pipette	Cotton wool stick	Screw cap

Contents of packagings

10 units:

- 10 pcs. Processing tube I, filled with 3,5 ml Medium A SAF solution
- 10 pcs. Processing tube II, with integrated filter-system
- 10 pcs. Sample spoon
- 10 pcs. Screw cap
- 10 pcs. One-way pipette 1 ml
- 10 pcs. Cotton wool stick
- 1 pc. 30 ml Combi Medium in brown glass
- 1 pc. Graduated dropper

50 units:









- 50 pcs. Processing tube I, filled with 3,5 ml Medium A SAF solution
- 50 pcs. Processing tube II, with integrated filter-system
- 50 pcs. Sample spoon
- 50 pcs. Screw cap
- 50 pcs. One-way pipette 1 ml
- 50 pcs. Cotton wool stick
- 1 pc. 75 ml Combi Medium in brown glass
- 1 pc. Graduated dropper

100 units:

- 100 pcs. Processing tube I, filled with 3,5 ml Medium A SAF solution
- 100 pcs. Processing tube II, with integrated filter-system
- 100 pcs. Sample spoon
- 100 pcs. Screw cap
- 100 pcs. One-way pipette 1 ml
- 100 pcs. Cotton wool stick
- 1 pc. 150 ml Combi Medium in brown glass
- 1 pc. Graduated dropper

SMALL ANIMALS

ParasiTrap® Premium Combipack Comfort System Bailenger

Order-No.: 7490 Bailenger	(10, 50, 100 pieces/box)	Premium Combi System	
			
Processing tube I (filled with 3,5 ml Medium A Bailenger solution)	Processing tube II (with integrated filter system)	Brown glass (with graduated dropper)	Combi Medium (incl. the sufficient quantity of media)
			
Sample spoon	One-way pipette	Cotton wool stick	Screw cap

Contents of packagings

10 units:

- 10 pcs. Processing tube I, filled with 3,5 ml Medium A Bailenger solution
- 10 pcs. Processing tube II, with integrated filter-system
- 10 pcs. Sample spoon
- 10 pcs. Screw cap
- 10 pcs. One-way pipette 1 ml
- 10 pcs. Cotton wool stick
- 1 pc. 30 ml Combi Medium in brown glass
- 1 pc. Graduated dropper

50 units:











- 50 pcs. Processing tube I, filled with 3,5 ml Medium A Bailenger solution
- 50 pcs. Processing tube II, with integrated filter-system
- 50 pcs. Sample spoon
- 50 pcs. Screw cap
- 50 pcs. One-way pipette 1 ml
- 50 pcs. Cotton wool stick
- 1 pc. 75 ml Combi Medium in brown glass
- 1 pc. Graduated dropper

100 units:

- 100 pcs. Processing tube I, filled with 3,5 ml Medium A Bailenger solution
- 100 pcs. Processing tube II, with integrated filter-system
- 100 pcs. Sample spoon
- 100 pcs. Screw cap
- 100 pcs. One-way pipette 1 ml
- 100 pcs. Cotton wool stick
- 1 pc. 150 ml Combi Medium in brown glass
- 1 pc. Graduated dropper

BIG ANIMALS

ParasiTrap® Flotation PLUS II BA DELUXE














Order-No.: 7491 FLOT II BA		(10 pieces/box)	„DELUXE“ Complete System	
				
Processing tube I (unfilled)		Processing tube II (with integrated filter system)		FixSepar® Flotation PLUS Medium (in brown glass with graduated dropper)
				
Shovel		Sample spoon		Screw cap
				
Plastic sample vessel (prefilled)	Funnel	Filter	Plastic vessel	

Contents of packagings

- 10 pcs. Processing tube I
- 10 pcs. Processing tube II, with integrated filter-system
- 10 pcs. Sample spoon
- 10 pcs. Screw cap
- 10 pcs. Shovel
- 10 pcs. Funnel
- 10 pcs. One-way filter
- 10 pcs. Plastic vessel
- 10 pcs. Plastic sample vessel with 30 ml FixSepar® Flotation PLUS Medium
- 1 pc. 30 ml FixSepar® Flotation PLUS Medium in brown glass with graduated dropper

BIG ANIMALS

ParasiTrap[®] ECO Big Animal DELUXE














Order-No.: 7492 ECO BA		(10 pieces/box)	„DELUXE“ Complete System	
				
Processing tube I (filled with 2 ml Medium A FixSepar [®] ECO solution)	Processing tube II (with integrated filter system)	Graduated dropper	Combi Medium	
				
Shovel	Sample spoon	Screw cap	2 one-way pipettes (1 and 3 ml)	
				
Plastic sample vessel (prefilled)	Funnel	Filter	Plastic vessel	
				
Cotton wool stick				

Contents of packagings

- 10 pcs. Processing tube I, filled with 2 ml Medium A FixSepar[®] ECO solution
- 10 pcs. Processing tube II, with integrated filter-system
- 10 pcs. Sample spoon
- 10 pcs. Screw cap
- 10 pcs. Shovel
- 10 pcs. Funnel
- 10 pcs. One-way pipette 1 ml
- 10 pcs. One-way pipette 3 ml
- 10 pcs. One-way filter
- 10 pcs. Plastic vessel
- 10 pcs. Plastic sample vessel with 30 ml Medium A FixSepar[®] ECO solution
- 10 pcs. Cotton wool stick
- 1 pc. 30 ml Combi Medium in brown glass
- 1 pc. Graduated dropper

BIG ANIMALS

ParasiTrap® AF Big Animal DELUXE














Order-No.: 7493 AF BA		(10 pieces/box)	„DELUXE“ Complete System	
				
Processing tube I (filled with 2 ml Medium AF solution)	Processing tube II (with integrated filter system)	Graduated dropper	Combi Medium	
				
Shovel	Sample spoon	Screw cap	2 one-way pipettes (1 and 3 ml)	
				
Plastic sample vessel (prefilled)	Funnel	Filter	Plastic vessel	
				
Cotton wool stick				

Contents of packaging

- 10 pcs. Processing tube I, filled with 2 ml Medium AF solution
- 10 pcs. Processing tube II, with integrated filter-system
- 10 pcs. Sample spoon
- 10 pcs. Screw cap
- 10 pcs. Shovel
- 10 pcs. Funnel
- 10 pcs. One-way pipette 1 ml
- 10 pcs. One-way pipette 3 ml
- 10 pcs. One-way filter
- 10 pcs. Plastic vessel
- 10 pcs. Plastic sample vessel with 30 ml Medium AF solution
- 10 pcs. Cotton wool stick
- 1 pc. 30 ml Combi Medium in brown glass
- 1 pc. Graduated dropper

BIG ANIMALS

ParasiTrap® SAF Big Animal DELUXE














Order-No.: 7494 SAF BA		(10 pieces/box)	„DELUXE“ Complete System	
				
Processing tube I (filled with 2 ml Medium A SAF solution)	Processing tube II (with integrated filter system)	Graduated dropper	Combi Medium	
				
Shovel	Sample spoon	Screw cap	2 one-way pipettes (1 and 3 ml)	
				
Plastic sample vessel (prefilled)	Funnel	Filter	Plastic vessel	
				
Cotton wool stick				

Contents of packaging

- 10 pcs. Processing tube I, filled with 2 ml Medium A SAF solution
- 10 pcs. Processing tube II, with integrated filter-system
- 10 pcs. Sample spoon
- 10 pcs. Screw cap
- 10 pcs. Shovel
- 10 pcs. Funnel
- 10 pcs. One-way pipette 1 ml
- 10 pcs. One-way pipette 3 ml
- 10 pcs. One-way filter
- 10 pcs. Plastic vessel
- 10 pcs. Plastic sample vessel with 30 ml Medium A SAF solution
- 10 pcs. Cotton wool stick
- 1 pc. 30 ml Combi Medium in brown glass
- 1 pc. Graduated dropper

BIG ANIMALS

ParasiTrap® Bailenger Big Animal DELUXE

Order-No.: 7495 Bailenger BA		(10 pieces/box)	„DELUXE“ Complete System	
				
Processing tube I (filled with 2 ml Medium A Bailenger solution)	Processing tube II (with integrated filter system)	Graduated dropper	Combi Medium	
				
Shovel	Sample spoon	Screw cap	2 one-way pipettes (1 and 3 ml)	
				
Plastic sample vessel (prefilled)	Funnel	Filter	Plastic vessel	
				
Cotton wool stick				

Contents of packaging

- 10 pcs. Processing tube I, filled with 2 ml Medium A Bailenger solution
- 10 pcs. Processing tube II, with integrated filter-system
- 10 pcs. Sample spoon
- 10 pcs. Screw cap
- 10 pcs. Shovel
- 10 pcs. Funnel
- 10 pcs. One-way pipette 1 ml
- 10 pcs. One-way pipette 3 ml
- 10 pcs. One-way filter
- 10 pcs. Plastic vessel
- 10 pcs. Plastic sample vessel with 30 ml Medium A Bailenger solution
- 10 pcs. Cotton wool stick
- 1 pc. 30 ml Combi Medium in brown glass
- 1 pc. Graduated dropper

ECONOMY VERSION FOR THE BIG ANIMAL DIAGNOSTICS

The Economy versions of 7491 FLOT II BA, 7492 ECO BA, 7493 AF BA, 7494 SAF BA, 7495 Bailenger BA contains only three shovels and funnels. The used shovels and funnels are re-usable after washing and disinfection. This alternative protects the environment and saves your budget.

If you need more information, don't hesitate, use our hotline:
+49 (0)8631 167 4855

THE BETTER TEAM WINS...



TOGETHER YOU HAVE THE BEST CHANCE







**The ParasiTrap[®] flotation- and sedimentation systems
are the best choice for your decision.**

 biosepar *HUMAN line*

**THE MODERN PARASITOLOGICAL DIAGNOSTICS
FOR THE HUMAN MEDICINE**



ParasiTrap® Flotation PLUS I HU

Order-No.: 8486 FLOT I HU		(10, 50, 100 pieces/box)	Complete system
			
Processing tube I	Processing tube II	FixSepar® Flotation PLUS Medium	FixSepar® Flotation PLUS Medium
(filled with 5,5 ml Medium A FixSepar® Flotation solution)	(with integrated filter system)	(in brown glass with graduated dropper)	(incl. the sufficient quantity of media)
			
Sample spoon		Screw cap	

Contents of packagings

10 units:

- 10 pcs. Processing tube I, filled with 5,5 ml FixSepar® Flotation PLUS Medium
- 10 pcs. Processing tube II, with integrated filter-system
- 10 pcs. Sample spoon
- 10 pcs. Screw cap
- 1 pc. 30 ml FixSepar® Flotation PLUS Medium in brown glass with graduated dropper









50 units:

- 50 pcs. Processing tube I, filled with 5,5 ml FixSepar® Flotation PLUS Medium
- 50 pcs. Processing tube II, with integrated filter-system
- 50 pcs. Sample spoon
- 50 pcs. Screw cap
- 1 pc. 125 ml FixSepar® Flotation PLUS Medium
- 1 pc. 30 ml FixSepar® Flotation PLUS Medium in brown glass with graduated dropper









100 units:

- 100 pcs. Processing tube I, filled with 5,5 ml FixSepar® Flotation PLUS Medium
- 100 pcs. Processing tube II, with integrated filter-system
- 100 pcs. Sample spoon
- 100 pcs. Screw cap
- 1 pc. 300 ml FixSepar® Flotation PLUS Medium
- 1 pc. 30 ml FixSepar® Flotation PLUS Medium in brown glass with graduated dropper

ParasiTrap[®] Human Premium Combipack Comfort System FixSepar[®] ECO









Order-No.: 8486 ECO	(50 pieces/box)*	Premium Combi System	
			
Processing tube I (filled with 3,5 ml Medium A FixSepar [®] ECO solution)	Processing tube II (with integrated filter system)	Brown glass (with graduated dropper)	Combi Medium (incl. the sufficient quantity of media)
			
Sample spoon	One-way pipette	Cotton wool stick	Screw cap

ParasiTrap[®] Human Premium Combipack Comfort System AF









Order-No.: 8486 AF	(50 pieces/box)*	Premium Combi System	
			
Processing tube I (filled with 3,5 ml Medium AF solution)	Processing tube II (with integrated filter system)	Brown glass (with graduated dropper)	Combi Medium (incl. the sufficient quantity of media)
			
Sample spoon	One-way pipette	Cotton wool stick	Screw cap

* By request: separated supply into 50 or 100 units

ParasiTrap[®] Human Premium CombiPack Comfort System SAF





Order-No.: 8486 SAF	(50 pieces/box)*	Premium Combi System	
			
Processing tube I (filled with 3,5 ml Medium A SAF solution)	Processing tube II (with integrated filter system)	Brown glass (with graduated dropper)	Combi Medium (incl. the sufficient quantity of media)
			
Sample spoon	One-way pipette	Cotton wool stick	Screw cap

ParasiTrap[®] Human Premium CombiPack Comfort System Bailenger





Order-No.: 8486 Bailenger	(50 pieces/box)*	Premium Combi System	
			
Processing tube I (filled with 3,5 ml Medium A Bailenger solution)	Processing tube II (with integrated filter system)	Brown glass (with graduated dropper)	Combi Medium (incl. the sufficient quantity of media)
			
Sample spoon	One-way pipette	Cotton wool stick	Screw cap

* By request: separated supply into 50 or 100 units





ParasiTrap[®] Super Combipack Comfort System FixSepar[®] ECO

Order-No.: 8484 ECO	(50 pieces/box)	Super Combi System	
			
Processing tube I (filled with 3,5 ml Medium A FixSepar [®] ECO solution)	Processing tube II (with integrated filter-system)	Brown glass (with graduated dropper)	Combi Medium (incl. the sufficient quantity of media)





ParasiTrap[®] Super Combipack Comfort System AF

Order-No.: 8484 AF	(50 pieces/box)	Super Combi System	
			
Processing tube I (filled with 3,5 ml Medium AF solution)	Processing tube II (with integrated filter-system)	Brown glass (with graduated dropper)	Combi Medium (incl. the sufficient quantity of media)








ParasiTrap® Super Combipack Comfort System SAF

Order-No.: 8484 SAF		(50 pieces/box)	Super Combi System	
				
Processing tube I (filled with 3,5 ml Medium A SAF solution)	Processing tube II (with integrated filter-system)	Brown glass (with graduated dropper)	Combi Medium (incl. the sufficient quantity of media)	








ParasiTrap® Super Combipack Comfort System Bailenger

Order-No.: 8484 Bailenger		(50 pieces/box)	Super Combi System	
				
Processing tube I (filled with 3,5 ml Medium A Bailenger solution)	Processing tube II (with integrated filter-system)	Brown glass (with graduated dropper)	Combi Medium (incl. the sufficient quantity of media)	

ParasiTrap® DUAL Diagnostic System for the parasitology



Order-No.: 8486 DUAL		(2x25 pieces/box)*		DUAL	
		+			
		+			
Processing tube I filled with		Processing tube II		Combi Medium	
25 pcs. Medium A Bailenger solution		25/25 (with integrated filter-system)		(incl. the sufficient quantity of media)	
					
Screw cap		One-way pipette		Cotton wool stick	

ParasiTrap® DUAL ECO Diagnostik System for the parasitology



Order-No.: 8486 DUAL ECO		(2x25 pieces/box)*		DUAL ECO	
		+			
		+			
Processing tube I filled with		Processing tube II		Combi Medium	
25 pcs. Medium A Bailenger solution		25/25 (with integrated filter-system)		(incl. the sufficient quantity of media)	
					
Screw cap		One-way pipette		Cotton wool stick	

* By request: separated supply into 50 or 100 units



ParasiTrap® FixSepar® ECO System

Order-No.: 8488 ECO	(50 pieces/box)	Complete system
		
Processing tube I (filled with 3,5 ml Medium A FixSepar® ECO)		Processing tube II (with integrated filter-system)



ParasiTrap® AF System

Order-No.: 8488 AF	(50 pieces/box)	Complete system
		
Processing tube I (filled with 3,5 ml Medium AF)		Processing tube II (with integrated filter-system)


ParasiTrap® SAF System

Order-No.: 8488 SAF	(50 pieces/box)	Complete system
		
Processing tube I (filled with 3,5 ml Medium A SAF)		Processing tube II (with integrated filter-system)


ParasiTrap® Bailenger System

Order-No.: 8488 Bailenger	(50 pieces/box)	Complete system
		
Processing tube I (filled with 3,5 ml Medium A Bailenger)		Processing tube II (with integrated filter-system)


ParasiTrap® Single tube FixSepar® ECO

Order-No.: 8490 ECO	(50 pieces/box)	Processing tube I
Processing tube I (filled with 3,5 ml Medium A FixSepar® ECO)		


ParasiTrap® Single tube AF

Order-No.: 8490 AF	(50 pieces/box)	Processing tube I
Processing tube I (filled with 3,5 ml Medium AF)		


ParasiTrap® Single tube SAF



Order-No.: 8490 SAF	(50 pieces/box)	Processing tube I
Processing tube I (filled with 3,5 ml Medium A SAF)		


ParasiTrap® Single tube Bailenger

Order-No.: 8490 Bailenger	(50 pieces/box)	Processing tube I
<p>Processing tube I (filled with 3,5 ml Medium A Bailenger)</p>		

ParasiTrap® Empty system


Order-No.: 8489	(50 pieces/box)	Processing tube I (single)
<p>Processing tube I (unfilled)</p>		

Order-No.: 8487	(50 pieces/box)	Unfilled system
		
<p>Processing tube I (unfilled)</p>	<p>Processing tube II (with integrated filter-system)</p>	


Order-No.: 8491	(50 pieces/box)	Processing tube II (single)
<p>Processing tube II (with integrated filter-system) Colours: black, green or blue</p>		

Biosepar solutions for the ParasiTrap[®] diagnostic system


Order-No.: 8487/AF, 8487/A SAF, 8487/A Bailenger	Medium A
Transport, processing and fixative solution	

	8487/AF	Transport, processing and fixative solution with formaldehyde
	8487/A SAF	Fixative and processing solution with formaldehyde
	8487/A Bailenger	Processing solution


Order-No.: 8487/B	Medium B
Separation solution	

	Ethyl Acetate
--	----------------------

Order-No.: 8487/C	Medium C
Staining solution	

	ParasiTrap Lugol (special modified and made for parasitological purposes)
---	---

Order-No.: 8487/CM	Combi Medium
Combined separation and staining medium	

	Plastic bottle filled with the new developed separation and staining solution for the parasitology. We recommend in additional the brown glass with graduated dropper (Order-No.: 8487/E 30).
---	---

Medium A FixSepar[®] ECO available only in the prefilled Processing tube I.

Solutions for the **ParasiTrap**[®] diagnostic system

Order-No.: 8500

FixSepar[®] Flotation PLUS Medium



Order-No.: 8535

Dispersion solution (for sample dilution)




Order-No.: 8536

Isopropanol 70% (for disinfection)





Order-No.: 8514


Hemolyse reagent









Dispenser

Order-No.: 8495/A	(1 piece/box)	Dispenser
<p>2,0-10,0 ml (for Medium AF/SAF)</p>		


Order-No.: 8495/B	(1 piece/box)	Dispenser
<p>1,0-5,0 ml (for Medium B)</p>		

Order-No.: 8495/C	(1 piece/box)	Dispenser
<p>0,1-1,0 ml (for Medium C)</p>		


Additional accessories for stool samples

Order-No.: 8492	(50 pieces/box)	Sample spoon
for stool samples human and small animals		
Order-No.: 8492 BA	(50 pieces/box)	Sample spoon
for stool samples big animals		
Order-No.: 8493/1	(50 pieces/box)	One-way pipette
for liquid stool samples collecting and processing in the ParasiTrap® system (1 ml)		
Order-No.: 8493/3	(50 pieces/box)	One-way pipette
for liquid stool samples collecting and processing in the ParasiTrap® system (3 ml)		
Order-No.: 8494	(50 pieces/box)	Cotton wool stick
for ParasiTrap® processing		
Order-No.: 8511	(10 pieces/box)	Funnel
additional accessory for 7491 FLOT II BA 7492 ECO BA 7493 AF BA 7494 SAF BA 7495 Bailenger BA		


Order-No.: 8512	(10 pieces/box)	Shovel
------------------------	-----------------	---------------

additional accessory for 7491 FLOT II BA 7492 ECO BA 7493 AF BA 7494 SAF BA 7495 Bailenger BA		
--	--	--


Order-No.: 8502	(1 piece/box)	Test tube rack
------------------------	---------------	-----------------------

for ParasiTrap [®] processing tube I and II (place for 25 tubes)		
---	--	--


Order-No.: 8503	(1 piece/box)	Test tube rack
------------------------	---------------	-----------------------


for 50 ml plastic vessel in the big animal diagnostics (24 holes) colours: yellow, blue, red, beige		
--	--	--

Order-No.: 8504	(1 piece/box)	Shakers vortex mixer
------------------------	---------------	-----------------------------


for Biosepar processing tubes		
-------------------------------	--	--

Additional accessories for the centrifugation


Order-No.: 8505	(1 piece/box)	Adapter for Biosepar tubes
Processing tube I and II adapter for centrifuge insert		


Order-No.: 8510	(1 piece/box)	Centrifuge insert
for Biosepar processing tubes		


Additional accessories for the flotation

Order-No.: 8501	(1 piece/box)	Wire loop
for the flotation method		


Additional accessories for the microscopy

Order-No.: 8506	(50 pieces/box)	Sample holder
with opaque field 0,8-1,0 mm thickness		




Order-No.: 8507	(50 pieces/box)	Covering glass
24 x 36 mm		

Order-No.: 8509	(50 pieces/box)	Covering glass
24 x 60 mm		

Screw cap for the Processing tube I and II

Order-No.: 8482	Screw cap
Additional screw cap for small centrifuge and quality control	

Brown glasses

Order-No.: 8487 E 30	(1 piece/box)	Dropper
<p>Brown glass (30 ml) with graduated dropper for the ParasiTrap® Combi Medium and Flotation PLUS Medium.</p>		
Order-No.: 8496/A	(1 piece/box)	Brown glass
<p>1000 ml for Medium A</p>		
Order-No.: 8496/B	(1 piece/box)	Brown glass
<p>500 ml for Medium B</p>		
Order-No.: 8496/C	(1 piece/box)	Brown glass
<p>250 ml for Medium C</p>	